

Product Datasheet

Human PXDN ELISA Kit (orb441851)

| | |
|--------------------------|---|
| Description | Human Peroxidasin Homolog (PXDN) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | Q92626 |
| Sensitivity | 0.059ng/mL |
| Entrez | 7837 |
| Expiration Date | Please enquire. |