

Product Datasheet

Human PYY3 ELISA Kit (orb441859)

Description	Human Peptide YY 3 (PYY3) competitive inhibition ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q5JQD4
Sensitivity	4.12pg/mL
Entrez	644059
Expiration Date	Please enquire.