

Product Datasheet

Canine RANTES ELISA Kit (orb441889)

Description	Canine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) Sandwich ELISA Kit
Reactivity	Canine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8HYS0
Sensitivity	7.8pg/mL
Entrez	403522
Expiration Date	Please enquire.