

Product Datasheet

Guinea pig RANTES ELISA Kit (orb441890)

Description	Guinea pig Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) Sandwich ELISA K...
Reactivity	Guinea pig
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P97272
Sensitivity	0.059ng/mL
Entrez	100135495
Expiration Date	Please enquire.