

Product Datasheet

Porcine RANTES ELISA Kit (orb441893)

Description	Porcine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) Sandwich ELISA Kit
Reactivity	Porcine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2EN89
Sensitivity	1.15ng/mL
Entrez	396613
Expiration Date	Please enquire.