

## Product Datasheet

### Human RIPK1 ELISA Kit (orb441953)

<b>Description</b>	Human Receptor Interacting Serine Threonine Kinase 1 (RIPK1) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q13546</a>
<b>Sensitivity</b>	0.057ng/mL
<b>Entrez</b>	<a href="#">8737</a>
<b>Expiration Date</b>	Please enquire.