

Product Datasheet

Mouse Reelin ELISA Kit (orb441957)

Description	Mouse Reelin (RL) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q60841
Sensitivity	6.5pg/mL
Entrez	19699
Expiration Date	Please enquire.