

Product Datasheet

Human RLN ELISA Kit (orb441960)

Description	Human Relaxin (RLN) competitive inhibition ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04808
Sensitivity	9.66pg/mL
Entrez	6013
Expiration Date	Please enquire.