

## Product Datasheet

### Rat RLN ELISA Kit (orb441963)

<b>Description</b>	Rat Relaxin (RLN) competitive inhibition ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.25-2000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<b>P01347</b>
<b>Sensitivity</b>	9.16pg/mL
<b>Entrez</b>	<b>25616</b>
<b>Expiration Date</b>	Please enquire.