

Product Datasheet

Mouse RSP01 ELISA Kit (orb442003)

Description	Mouse R-Spondin 1 (RSP01) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z132
Sensitivity	5.8pg/mL
Entrez	192199
Expiration Date	Please enquire.