

Product Datasheet

Human S100A12 ELISA Kit (orb442024)

Description	Human S100 Calcium Binding Protein A12 (S100A12) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P80511
Sensitivity	34pg/mL
Entrez	6283
Expiration Date	Please enquire.