

Product Datasheet

Canine SCF ELISA Kit (orb442079)

Description	Canine Stem Cell Factor (SCF) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q06220
Sensitivity	6.5pg/mL
Entrez	403507
Expiration Date	Please enquire.