

Product Datasheet

Human SCNN1a ELISA Kit (orb442094)

| | |
|--------------------------|--|
| Description | Human Amiloride Sensitive Sodium Channel Subunit Alpha (SCNN1a) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | P37088 |
| Sensitivity | 0.067ng/mL |
| Entrez | 6337 |
| Expiration Date | Please enquire. |