

Product Datasheet

General SCT ELISA Kit (orb442098)

| | |
|--------------------------|--|
| Description | General Salmon Calcitonin (SCT) competitive inhibition ELISA Kit |
| Reactivity | All |
| Range | 0.312-20ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | biological agents. |
| Sample Types | biological agents. |
| Sensitivity | 0.085ng/mL |
| Expiration Date | Please enquire. |