

Product Datasheet

Human SECTM1 ELISA Kit (orb442120)

Description	Human Secreted And Transmembrane Protein 1 (SECTM1) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, urine and other biological fluids.
Sample Types	tissue homogenates, urine and other biological fluids.
Uniprot ID	Q8WVN6
Sensitivity	0.27ng/mL
Entrez	6398
Expiration Date	Please enquire.