

Product Datasheet

Canine SELE ELISA Kit (orb442122)

Description	Canine Selectin, Endothelium (SELE) Sandwich ELISA Kit
Reactivity	Canine
Range	6.25-400ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P33730
Sensitivity	2.38ng/mL
Entrez	403999
Expiration Date	Please enquire.