

Product Datasheet

Mouse SHBG ELISA Kit (orb442193)

Description	Mouse Sex Hormone Binding Globulin (SHBG) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P97497
Sensitivity	12.6pg/mL
Entrez	20415
Expiration Date	Please enquire.