

Product Datasheet

Rat SHBG ELISA Kit (orb442194)

Description	Rat Sex Hormone Binding Globulin (SHBG) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P08689
Sensitivity	6.7pg/mL
Entrez	24775
Expiration Date	Please enquire.