

Product Datasheet

Human SIGLEC10 ELISA Kit (orb442206)

Description	Human Sialic Acid Binding Ig Like Lectin 10 (SIGLEC10) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q96LC7
Sensitivity	0.062ng/mL
Entrez	89790
Expiration Date	Please enquire.