

Product Datasheet

Human SIGLEC8 ELISA Kit (orb442209)

Description	Human Sialic Acid Binding Ig Like Lectin 8 (SIGLEC8) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	cell lysates.
Sample Types	cell lysates.
Uniprot ID	Q9NYZ4
Sensitivity	0.058ng/mL
Entrez	27181
Expiration Date	Please enquire.