

Product Datasheet

Mouse SIRT6 ELISA Kit (orb442224)

Description	Mouse Sirtuin 6 (SIRT6) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P59941
Sensitivity	0.063ng/mL
Entrez	50721
Expiration Date	Please enquire.