



## Product Datasheet

### CD46 antibody (orb44229)

<b>Description</b>	Mouse Monoclonal to CD46.
<b>Reactivity</b>	Bovine, Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IP, WB
<b>Immunogen</b>	HPB-ALL human T cell line
<b>Target</b>	CD46
<b>Preservatives</b>	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 2-8°C. Do not freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	Flow cytometry: Recommended dilution: 1 µg/ml. Western blotting: Non-reducing conditions.
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified by protein-A affinity chromatography.
<b>Uniprot ID</b>	<b>P15529</b>
<b>Entrez</b>	<b>4179</b>
<b>Dilution Range</b>	Flow cytometry: Recommended dilution: 1 µg/ml. Western blotting: Non-reducing conditions.
<b>Expiration Date</b>	12 months from date of receipt.



Separation of human lymphocytes (red-fil...



Flow cytometry surface staining pattern ...