

Product Datasheet

Human SOSTDC1 ELISA Kit (orb442305)

Description	Human Sclerostin Domain Containing Protein 1 (SOSTDC1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q6X4U4
Sensitivity	0.059ng/mL
Entrez	25928
Expiration Date	Please enquire.