

Product Datasheet

Human SPB ELISA Kit (orb442322)

Description	Human Surfactant Associated Protein B (SPB) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, lung lavage fluid, cell culture supernatants and other biological fluids.
Uniprot ID	P07988
Sensitivity	0.54ng/mL
Entrez	6439
Expiration Date	Please enquire.