

Product Datasheet

Mouse SPD ELISA Kit (orb442328)

Description	Mouse Surfactant Associated Protein D (SPD) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.
Uniprot ID	P50404
Sensitivity	0.285ng/mL
Entrez	20390
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)