

Product Datasheet

Human STS ELISA Kit (orb442402)

| | |
|--------------------------|--|
| Description | Human Steroid Sulfatase Isozyme S (STS) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 4.687-300ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | tissue homogenates, cell lysates and other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | P08842 |
| Sensitivity | 1.9ng/mL |
| Entrez | 412 |
| Expiration Date | Please enquire. |