

Product Datasheet

Human TBG ELISA Kit (orb442462)

Description	Human Thyroxine Binding Globulin (TBG) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05543
Sensitivity	0.33ng/mL
Entrez	6906
Expiration Date	Please enquire.