

## Product Datasheet

### Human TERF2IP ELISA Kit (orb442481)

<b>Description</b>	Human Telomeric Repeat Binding Factor 2 Interacting Protein (TERF2IP) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates and other biological fluids.
<b>Sample Types</b>	tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9NYB0</a>
<b>Sensitivity</b>	0.66ng/mL
<b>Entrez</b>	<a href="#">54386</a>
<b>Expiration Date</b>	Please enquire.