

Product Datasheet

Human TFF1 ELISA Kit (orb442487)

| | |
|--------------------------|--|
| Description | Human Trefoil Factor 1 (TFF1) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 125-8000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | P04155 |
| Sensitivity | 54pg/mL |
| Entrez | 7031 |
| Expiration Date | Please enquire. |