

## Product Datasheet

### Mouse TFF2 ELISA Kit (orb442490)

<b>Description</b>	Mouse Trefoil Factor 2 (TFF2) Sandwich ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25-2000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9QX97</a>
<b>Sensitivity</b>	12.5pg/mL
<b>Entrez</b>	<a href="#">21785</a>
<b>Expiration Date</b>	Please enquire.