

Product Datasheet

Porcine Tissue Factor ELISA Kit (orb442497)

Description	Porcine Tissue Factor (TF) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q6RCS8
Sensitivity	0.061ng/mL
Entrez	396677
Expiration Date	Please enquire.