

Product Datasheet

Human TFPI2 ELISA Kit (orb442498)

Description	Human Tissue Factor Pathway Inhibitor 2 (TFPI2) Sandwich ELISA Kit
Reactivity	Human
Range	6.25-400pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P48307
Sensitivity	2.75pg/mL
Entrez	7980
Expiration Date	Please enquire.