

Product Datasheet

Rat TFPI ELISA Kit (orb442501)

Description	Rat Tissue Factor Pathway Inhibitor (TFPI) Sandwich ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma, urine, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	plasma, urine, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	Q02445
Sensitivity	0.66ng/mL
Entrez	29436
Expiration Date	Please enquire.