

Product Datasheet

Human Tie1 ELISA Kit (orb442549)

Description	Human Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) Sandwich ELISA ...
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P35590
Sensitivity	0.62ng/mL
Entrez	7075
Expiration Date	Please enquire.