

Product Datasheet

Mouse Tie1 ELISA Kit (orb442550)

Description	Mouse Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) Sandwich ELISA ...
Reactivity	Mouse
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q06806
Sensitivity	0.247ng/mL
Entrez	21846
Expiration Date	Please enquire.