

Product Datasheet

Human TKTL1 ELISA Kit (orb442597)

| | |
|--------------------------|------------------------------------------------------------------------------------------|
| Description | Human Transketolase Like Protein 1 (TKTL1) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 0.156-10ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | tissue homogenates, cell lysates and other biological fluids. |
| Sample Types | tissue homogenates, cell lysates and other biological fluids. |
| Uniprot ID | P51854 |
| Sensitivity | 0.062ng/mL |
| Entrez | 8277 |
| Expiration Date | Please enquire. |