

Product Datasheet

Rat TLR2 ELISA Kit (orb442608)

Description	Rat Toll Like Receptor 2 (TLR2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	E9PTD9
Sensitivity	0.115ng/mL
Entrez	310553
Expiration Date	Please enquire.