

Product Datasheet

Mouse TLR6 ELISA Kit (orb442617)

| | |
|--------------------------|------------------------------------------------------------------------------------------|
| Description | Mouse Toll Like Receptor 6 (TLR6) Sandwich ELISA Kit |
| Reactivity | Mouse |
| Range | 0.312-20ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Sample Types | tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | Q9EPW9 |
| Sensitivity | 0.125ng/mL |
| Entrez | 21899 |
| Expiration Date | Please enquire. |