

## Product Datasheet

### Canine Thrombomodulin ELISA Kit (orb442622)

<b>Description</b>	Canine Thrombomodulin (TM) Sandwich ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	15.625-1000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	plasma.
<b>Sample Types</b>	plasma.
<b>Uniprot ID</b>	<a href="#">Q5W7P8</a>
<b>Sensitivity</b>	5.8pg/mL
<b>Entrez</b>	<a href="#">474355</a>
<b>Expiration Date</b>	Please enquire.