

Product Datasheet

Porcine Thrombomodulin ELISA Kit (orb442630)

Description	Porcine Thrombomodulin (TM) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	B3STX8
Sensitivity	0.29ng/mL
Entrez	100157642
Expiration Date	Please enquire.