

Product Datasheet

Human TMSb10 ELISA Kit (orb442637)

Description	Human Thymosin Beta 10 (TMSb10) competitive inhibition ELISA Kit
Reactivity	Human
Range	7.812-500ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P63313
Sensitivity	3.52ng/mL
Entrez	9168
Expiration Date	Please enquire.