

Product Datasheet

Porcine tPA ELISA Kit (orb442684)

Description	Porcine Plasminogen Activator, Tissue (tPA) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma and other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	Q85Q23
Sensitivity	0.32ng/mL
Entrez	397121
Expiration Date	Please enquire.