

Product Datasheet

Canine TPO ELISA Kit (orb442702)

Description	Canine Thrombopoietin (TPO) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P42705
Sensitivity	6.9pg/mL
Entrez	100855822
Expiration Date	Please enquire.