

Product Datasheet

Human TREM2 ELISA Kit (orb442730)

Description	Human Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9NZC2
Sensitivity	21.9pg/mL
Entrez	54209
Expiration Date	Please enquire.