

Product Datasheet

Human TRH ELISA Kit (orb442738)

Description	Human Thyrotropin Releasing Hormone (TRH) competitive inhibition ELISA Kit
Reactivity	Human
Range	7.812-500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P20396
Sensitivity	3.25pg/mL
Entrez	7200
Expiration Date	Please enquire.