

Product Datasheet

Mouse TRH ELISA Kit (orb442739)

Description	Mouse Thyrotropin Releasing Hormone (TRH) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	7.812-500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q62361
Sensitivity	2.77pg/mL
Entrez	22044
Expiration Date	Please enquire.