

Product Datasheet

Mouse TRPV2 ELISA Kit (orb442751)

Description	Mouse Transient Receptor Potential Cation Channel Subfamily V, Member 2 (TRPV2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9WTR1
Sensitivity	0.049ng/mL
Entrez	22368
Expiration Date	Please enquire.