

Product Datasheet

Canine TSH ELISA Kit (orb442759)

| | |
|--------------------------|--|
| Description | Canine Thyroid Stimulating Hormone (TSH) competitive inhibition ELISA Kit |
| Reactivity | Canine |
| Range | 46.875-3000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Sensitivity | 11.25pg/mL |
| Expiration Date | Please enquire. |