

Product Datasheet

General TXB2 ELISA Kit (orb442785)

Description	General Thromboxane B2 (TXB2) competitive inhibition ELISA Kit
Reactivity	All
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	7.13pg/mL
Expiration Date	Please enquire.