

Product Datasheet

Human UBE2I ELISA Kit (orb442801)

Description	Human Ubiquitin Conjugating Enzyme E2I (UBE2I) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P63279
Sensitivity	0.053ng/mL
Entrez	7329
Expiration Date	Please enquire.